

Patent

Docket No. 1759.15103

**EXPANDABLE PREFORMED STRUCTURES FOR
DEPLOYMENT IN INTERIOR BODY REGIONS**

FIELD OF THE INVENTION

5 The invention relates to expandable structures,
which, in use, are deployed in interior body regions of
humans and other animals.

BACKGROUND OF THE INVENTION

10 The deployment of expandable structures,
generically called "balloons," into cancellous bone is
known. For example, U.S. Patents 4,969,888 and 5,108,404
disclose apparatus and methods using expandable structures
in cancellous bone for the fixation of fractures or other
osteoporotic and non-osteoporotic conditions of human and
animal bones.

15 SUMMARY OF THE INVENTION

When deployed in cancellous bone, expandable
structures should undergo expansion and distention without
failure. Furthermore, such structures, when distended,
should generally match the geometry of the interior bone
space in which the structure is deployed. In addition, such
structures should allow preferential expansion to areas of
lowest bone density. Exposure to cancellous bone also
requires materials that exhibit superior resistance to
surface abrasion and tensile stresses.

25 It is has been discovered that expandable
structures made from an elastomer material, e.g.,
polyurethane, which have been preformed to a desired shape,
e.g., by exposure to heat and pressure, can undergo
controlled expansion and further distention in cancellous
bone, without failure, while exhibiting superior resistance

30

09037350-041501

DND
721

to surface abrasion and puncture when contacting cancellous bone.

Features and advantages of the inventions are set forth in the following Description and Drawings, as well as in the appended Claims.

BRIEF DESCRIPTION OF THE DRAWINGS

Fig. 1 is a coronal view of a vertebral body;

Fig. 2 is a lateral view of the vertebral body shown in Fig. 1;

Fig. 3 is a plan view of a tool which carries at its distal end an expandable structure that embodies features of the invention;

Figs. 4A and 4B are enlarged side views of the expandable structure carried by the tool shown in Fig. 3;

Fig. 5 is a perspective end view of a tube made of a polyurethane or elastomer material prior to being formed into the expandable structure shown in Fig. 4A;

Fig. 6 is a top perspective view of the tube shown in Fig. 5 positioned in a shape-forming fixture, of which parts are broken away to permit viewing its interior;

Fig. 7 is a top perspective view of the shape-forming fixture shown in Fig. 6, in use applying heat and pressure to a region of the tube to form a shaped, expandable region;

Fig. 8 is a coronal view of the vertebral body shown in Fig. 1, with the tool shown in Fig. 3 deployed to compress cancellous bone as a result of inflating the expandable structure;

Fig. 9 is a coronal view of the vertebral body shown in Fig. 8, upon removal of the tool, showing the cavity formed by the compression of cancellous bone by the expandable structure;

Fig. 10 is a graph which plots the effects of increasing pressure applied to the interior of the structure to the expanded volume of the structure;

Fig. 12 is a side view of a complex structure which includes several expandable segments spaced along its length;

Fig 14 is a top view of a kit which holds the tool shown in Fig. 3 in a sealed, sterile environment prior to use;

Fig. 16 is a side view, partly in section, of a composite expandable structure.

DETAILED DESCRIPTION OF THE PREFERRED EMBODIMENTS

The new systems and methods will be more
35 specifically described in the context of the treatment of

human vertebra. Of course, other human or animal bone types can be treated in the same or equivalent fashion.

Fig. 1 shows a coronal (top) view of a human lumbar vertebra 12. Fig. 2 shows a lateral (side) view of the vertebra. The vertebra 12 includes a vertebral body 26, which extends on the anterior (i.e., front or chest) side of the vertebra 12. The vertebral body 26 is shaped generally like a marshmallow.

As Figs. 1 and 2 show, the vertebral body 26 includes an exterior formed from compact cortical bone 28. The cortical bone 28 encloses an interior volume of reticulated cancellous, or spongy, bone 32 (also called medullary bone or trabecular bone).

The spinal canal 36 (see Fig. 1), is located on the posterior (i.e., back) side of each vertebra 12. The spinal cord (not shown) passes through the spinal canal 36. The vertebral arch 40 surrounds the spinal canal 36. Left and right pedicles 42 of the vertebral arch 40 adjoin the vertebral body 26. The spinous process 44 extends from the posterior of the vertebral arch 40, as do the left and right transverse processes 46.

It may be indicated, due to disease or trauma, to compress cancellous bone within the vertebral body. The compression, for example, can be used to form an interior cavity, which receives a filling material, e.g., a flowable material that sets to a hardened condition, like bone cement, allograft tissue, autograft tissue, hydroxyapatite, or synthetic bone substitute, as well as a medication, or combinations thereof, to provide improved interior support for cortical bone or other therapeutic functions, or both. The compaction of cancellous bone also exerts interior force upon cortical bone, making it possible to elevate or push broken and compressed bone back to or near its original prefracture, or other desired, condition.

I. Preformed Expandable Structures

Fig. 3 shows a tool 48 for accessing bone for the purpose of compacting cancellous bone. The tool 48 includes a catheter tube 50 having a proximal end 52 and a distal end 54. The proximal end 52 carries a handle 14 to facilitate gripping and maneuvering the tube 50. The proximal end 52 also carries a fitting 122 to enable connection of the tool 48 to external equipment, as will be described later. The distal end 54 of the tool 48 carries a structure 56, which, in use, is intended to be expanded in cancellous bone, e.g., in the vertebral body 26 shown in Figs. 1 and 2.

A. Desired Physical and Mechanical Properties

The material from which the structure 56 is made should possess various physical and mechanical properties to optimize its functional capabilities to compact cancellous bone. The three most important properties are the ability to expand its volume; the ability to deform in a desired way when expanding and assume a desired shape inside bone; and the ability to withstand abrasion, tearing, and puncture when in contact with cancellous bone.

1. Expansion Property

A first desired property for the structure material is the ability to expand or otherwise increase its volume without failure. This property enables the structure 56 to be deployed in a collapsed, low profile condition subcutaneously, e.g., through a cannula, into the targeted bone region. This property also enables the expansion of the structure 56 inside the targeted bone region to press against and compress surrounding cancellous bone, or move cortical bone to a prefracture or other desired condition, or both.

The expansion property for the material can be characterized, e.g., by ultimate elongation properties, which indicate the greatest degree of expansion that the material can accommodate prior to failure. An ultimate elongation of at least about 300% before material failure

5

2. Shape Property

10

15

35

achieve the desired therapeutic result. Such movement is not *per se* harmful, as that term is used in this Specification, because it is indicated to achieve the desired therapeutic result. By definition, harm results when expansion of the structure 56 results in a worsening of the overall condition of the bone and surrounding anatomic structures, for example, by injury to surrounding tissue or causing a permanent adverse change in bone biomechanics.

As one general consideration, in cases where the bone disease causing fracture (or the risk of fracture) is the loss of cancellous bone mass (as in osteoporosis), the selection of the expanded shape of the structure 56 inside bone should take into account the cancellous bone volume which should be compacted to achieve the desired therapeutic result. An exemplary range is about 30% to 90% of the cancellous bone volume, but the range can vary depending upon the targeted bone region. Generally speaking, compacting less of the cancellous bone volume leaves more uncompact, diseased cancellous bone at the treatment site.

Another general guideline for the selection of the expanded shape of the structure 56 inside bone is the amount that the targeted fractured bone region has been displaced or depressed. The controlled deformation diameter expansion of the structure 56 within the cancellous bone region inside a bone can elevate or push the fractured cortical wall back to or near its anatomic position occupied before fracture occurred. Generally speaking, inadequate compaction of cancellous bone results in less lifting of contiguous cortical bone.

For practical reasons, it is desired that the expanded shape of the structure 56 inside bone, when in contact with cancellous bone, substantially conforms to the shape of the structure 56 outside bone, when in an open air environment. This allows the physician to select in an open air environment a structure having an expanded shape desired to

meet the targeted therapeutic result, with the confidence that the expanded shape inside bone will be similar in important respects.

5 An optimal degree of shaping can be achieved by material selection and by special manufacturing techniques, e.g., thermoforming or blow molding, as will be described in greater detail later.

3. Toughness Property

10 A third desired property for the material of the structure 56 is the ability to resist surface abrasion, tearing, and puncture when in contact with cancellous bone.

15 This property can be characterized in various ways. For example, a Taber Abrasion Resistance Value of less than about 90 mg loss indicates resistance to puncture when contacting cancellous bone. A Rotating Drum Abrasion Resistance Value of less than 70 mm³ also indicates resistance to puncture when contacting cancellous bone. This property can further be characterized, e.g., by an Elmendorf tear strength of greater than about 280 lbf/in, which
20 indicates resistance to failure caused by cancellous bone abrasion. This property can also be characterized, e.g., by a Shore Hardness value of less than about 100 A. This value indicates a degree of elasticity, flexibility, and ductility.

25 Materials with a Taber Abrasion Resistance Value greater than about 90 mg loss, or a Rotating Drum Abrasion Resistance Value greater than about 70 mm³, or an Elmendorf tear strength value of less than about 280 lbf/in, or a Shore Hardness value greater than about 100 A are not well
30 suited for expansion in cancellous bone, because failure may occur prior to expansion to the desired diameter.

B. Enhanced Expansion and Shape Properties

35 The expansion and shape properties just described can be enhanced and further optimized for compacting cancellous bone by selecting an elastomer material, which also possess

09837350-041904

5

10

15

20

35

How it is formed

5

10

15

20

25

30

35

be circulated. The speed at which cooling occurs affects the time of the overall process.

After cooling, the application of pressurized fluid is terminated. The now preformed structure 56 is removed from the cavity 24.

The normally expanded shape characteristics of the structure 56 can be achieved by other techniques. For example, the structure 56 can be formed by dipping, lost wax casting, or injection molding.

Upon removal from the fixture 10, the structure 56 is secured to the distal end 54 of the catheter tube 50. The structure of the catheter tube 50 can vary and is not critical to the invention per se. The materials for the catheter tube 50 are selected to facilitate advancement of the structure 56 into an interior body region. The catheter tube 50 can be constructed, for example, using standard flexible, medical grade plastic materials, like vinyl, nylon, polyethylenes, ionomer, polyurethane, and polyethylene terephthalate (PET). The catheter tube 50 can also include more rigid materials to impart greater stiffness and thereby aid in its manipulation. More rigid materials that can be used for this purpose include Kevlar™ material, PEBAX™ material, stainless steel, nickel-titanium alloys (Nitinol™ material), and other metal alloys.

In the illustrated embodiment (as best shown in Fig. 4A), the catheter tube 50 includes an interior bore 60, in which an auxiliary tube 58 is secured. It should be appreciated that the catheter tube 50 can have more than a single interior lumen, and can, e.g., have an array of multiple lumens. In the illustrated embodiment, The auxiliary tube 58 extends through the interior bore 60 and beyond the distal end 54 of the catheter tube 50. One end region 18 of the tube 16 is secured to the distal end 54 of the catheter tube 50, while the other end region 20 is secured to the free extended end 62 of the auxiliary tube

58. The end regions 18 and 20 can be secured, e.g., using adhesive or thermal bonding processes.

By drawing a vacuum (i.e., negative pressure) inside the structure 56, resident air volume is removed, and the diameter of the region 22 is diminished from its normally expanded shape D3 to a substantially collapsed, and not inflated diameter D2. The collapsed diameter D2 is, due to forming during the heat and pressure shaping process, still different than the extruded or molded diameter D1. When substantially collapsed or not inflated, the structure 56 exhibits a low profile, ideal for insertion into the targeted cancellous bone region. The low profile can be further reduced to aid insertion, if desired, by enclosing the structure 56 within a constricted introducing sleeve, or by coating the structure 56 with a lubricious material, such as silicone, or both.

As Figs. 3 and 4 show, the interior bore 60 of the catheter tube 50 can be coupled (via the fitting 122) to a source 68 of fluid, for example, sterile saline, or a radiopaque contrast medium, which permits x-ray visualization of the structure 56. The interior bore 60 conveys the fluid into the region 22. The increase of volume within the region up to a given threshold amount (designated V(D3) in Fig. 10) will return the intermediate region 22 from the collapsed diameter D2 to the normal (i.e., enlarged, but not distended) geometry, having the shape and diameter D3.

~~When in its normally enlarged shape D3, the material of the structure 56 in the region 22 is not significantly stretched or stressed, because it has been previously expanded in a stress-relieved condition into this geometry in the cavity 24.~~

~~The~~ The magnitude of the radius of expansion (and thus diameter D3) depends upon the relative increase in diameter in the region 22 brought about by exposure to heat and

0522E550 041804

5 As Fig. 4B shows, due to expansion of heat-softened
material under pressure in the cavity 24, the wall thickness
of the structure 56 is not uniform. The region 22 has a
minimum wall thickness T3 when in its normally enlarged
diameter D3, which is less than the normal extruded or
10 molded wall thickness (T1) of the tube 86.

Continued volume flow of pressurized fluid into the structure 56 at the threshold pressure $P(t)$ continues to increase the interior volume of the structure 56. As its volume increases, the shaped region 22 of the structure 56 continues to enlarge beyond the normal diameter $D3$ toward a distended shape and geometry, designated $D4$ in Fig. 4. The threshold pressure $P(t)$ stays generally constant as volume increases between $D3$ and $D4$. As long as volume is controlled (i.e., so as not to substantially exceed $D4$), there is no need for an external pressure regulator. Volume expansion between $D3$ and $D4$ at a substantially constant pressure occurs because of the material properties of the structure 56, and not because of some external pressure control mechanism.

25 Enlargement of the structure in the region between D3
and D4 stretches the material in the region 22 beyond its
stress-relieved condition. Consequently, the wall thickness
T4 at the distended geometry D4 is less than the minimum
wall thickness T3 of the normally enlarged diameter D3.
30 However, the distended geometry generally maintains the
preformed shape dictated by the cavity 24 (which, in the
illustrated embodiment, is spherical).

In the expansion region between D3 and D4, the addition of fluid volume at substantially constant $P(t)$ stretches the material, causing the radius of the structure to increase

5

10

15

20

25

30

35

80 can be attached to the distal end of a catheter tube 50 in the same fashion as the structure 56 shown in Figs. 3 and 4.

5 The structure 80 possesses, in an open air environment, a normal expanded shape, having diameter D7 (shown in phantom lines in Fig.12). The normal shape and diameter D7 for the regions 82 and 84 generally correspond with the shape and dimension of the cavities 92 and 94, respectively.

10 When an interior vacuum is drawn, removing air from the structure 80, the structure 80 assumes a substantially collapsed, and not inflated geometry, shown in phantom lines D6 in Fig. 12. Due to the application of heat and pressure upon the region 88, diameter D6 for each region 82 and 84 is larger than the normal extruded or molded outside diameter
15 D5 of the original extruded tube 86.

20 The regions 82 and 84 are separated by a tubular neck 98, which segments the structure 80 into two expandable regions 82 and 84. When substantially collapsed under vacuum or not inflated, the structure 80 exhibits a low profile, ideal for insertion into the targeted cancellous bone region.

25 The introduction of fluid volume back into the tube 86 will cause each region 82 and 84 to return from the collapsed diameter D6 back to the normal, enlarged, but not distended geometry, having the shape and diameter shown in phantom lines D7 in Fig. 12.

30 In the illustrated embodiment, the first and second shaped regions 82 and 84 have generally the same radius of expansion and thus the same non-distended shape and diameter D7. Alternatively, each region 82 and 84 can have a different radius of expansion, and thus a different non-distended shape and diameter. Regardless, when in the normal, non-distended diameter D7, the material of the structure 80 in the region 88 is not significantly stretched
35 or stressed, because the regions 82 and 84 have been

TOP SECRET

5

10

15

20

35

provides multiple expandable regions 82 and 84 of the same or different enlarged outside diameters (D7), segmented by a neck region 98, in which $D6 > D5$; $D7 > D6$; and $D7 > D9$.

By continuing to apply fluid volume at a constant pressure at a threshold amount $P(t)$, and thereby increasing the volume within the structure 80, the shaped regions 82 and 84 of the structure 80 will continue to enlarge beyond diameter D7 to a distended shape and geometry, designated D8 in Fig. 12. The wall thickness T7 further decreases and approaches T8. As the regions 82 and 84 approach diameter D8, the diameter D9 of the neck region 98 will likewise increase toward diameter D10, as Fig. 12 shows, providing more uniform, elongated surface contact with cancellous bone.

Enlargement of the structure 80 beyond diameter D7 stretches the material in the regions 82, 84, and 98 beyond their stress-relieved condition, although the distended geometry of the regions 82 and 84 will, in important respects, maintain the preformed shape dictated by the cavities 92 and 94. As before explained in conjunction with the structure 56, the material in the regions 82 and 84 has already been stress-relieved in the desired shape at the normal diameter D7. As previously explained, enlargement toward the distended diameter D8 occurs at substantially constant pressure (as Fig. 10 exemplifies), and at increasing material stress.

The degree of stretching at a substantially constant incremental pressure condition can be tailored to achieve a desired, fully distended diameter D8. The final, fully distended diameter D8 can be selected to match the dimensions of the targeted cancellous bone region. The controlled stretching of the segmented regions 82 and 84 in tandem can provide an equal volume compression of cancellous bone with a major diameter that is less than a single non-segmented region (i.e., one without the neck region 98).

Stated another way, segmented regions 82 and 84, when expanded to a given inflation volume, have a diameter less than a sphere expanded to an equal inflation volume.

While expanding in the region between D7 and D8, the structure 80, like the structure 56, when inside bone, assumes an increasingly larger surface and volume, thereby compacting surrounding cancellous bone. Inflation in cancellous bone may occur at the same threshold pressure $P(t)$ as outside bone. However, an increase in the threshold inflation pressure $P(t)$ inside bone may be required, due to the density of the cancellous bone and resistance of the cancellous bone to compaction.

3. Composite Expandable Structures

In the previous embodiments, the material of the structure 56 or 80 is selected to integrate all desired physical and mechanical requirements of expansion, shape, and toughness. Fig. 16 exemplifies a composite expandable structure 130, in which the desired physical and mechanical requirements are segregated by the use of different materials.

As shown in Fig. 16, the composite structure 130 includes an inner expandable body 132 made of a first material that meets one or more of the desired requirements of expansion and shape. The composite structure 130 includes an outer expandable body or shell 134, which is made of a second material that meets the desired requirement of toughness. The shell 134 encapsulates and protects the inner expandable body 132 from surface abrasion, tearing, or puncture due to contact with cancellous bone.

The shell 134 can comprise a material applied to the surface of the inner body by various dipping, painting, or coating techniques. Alternatively, the shell 134 can comprise a bag or sock, into which the inner body 132 is placed prior to deployment. The material for the shell 134 can comprise, e.g., rubber, silicone, ethylene vinyl

0092750 04190
FOGTHO 092750

acetate, polyurethane, polyethylene, or multi-filament woven material or fabric or other polymer material.

The composite structure 130 makes it possible to isolate the expansion and shape requirements from the toughness requirement. A material completely or partially failing to meet the toughness requirement can nevertheless be used for the inner body 132 to optimize the expansion and shape requirements of the structure 130. The inner body 132 imparts its optimized expansion and shape characteristics to cancellous bone, while the shell 134 imparts its optimized toughness characteristic to the overall composite structure 130.

II. Deployment of Preformed Expandable Structures in Bone

The structure 56 or 80 or 130 can be inserted into bone in accordance with the teachings of U.S. Patents 4,969,888 and 5,108,404, which are incorporated herein by reference. As Fig. 8 shows, access can be accomplished, for example, by drilling an access portal 64 through a side of the vertebral body 26. This is called a postero-lateral approach. Alternatively, the access portal can pass through either pedicle 42, which called a transpedicular approach.

A guide sheath or cannula 66 is placed into communication with the access portal 64. The catheter tube 50 is advanced through the cannula 66 to deploy the structure (Fig. 8 shows structure 56) into contact with cancellous bone 32. The structure 56 is in its normally collapsed and not inflated condition (shown as phantom line diameter D2 in Fig. 8) during deployment. Access in this fashion can be accomplished using a closed, minimally invasive procedure or with an open procedure.

The materials for the catheter tube 50 are selected to facilitate advancement of the expandable structure 56 into cancellous bone 32. The catheter tube 50 can be constructed, for example, using standard flexible, medical grade plastic materials, like vinyl, nylon, polyethylenes,

5

10

20

25

35

the structure 56. The stylet 74 can be made, e.g., from stainless steel or molded plastic material. The presence of the stylet 74 serves to keep the structure 56 in the desired distally straightened condition during passage through the guide sheath 66 into the targeted bone region, as Fig. 8 shows.

As further shown in Fig. 11, the stylet 74 can have a preformed memory, to normally bend its distal region. The memory is overcome to straighten the stylet 14 when confined within the guide sheath 66. However, as the structure 56 and preformed stylet 74 advance free of the guide sheath 66 and pass into the targeted region, the preformed memory bends the stylet 74. The bending stylet 74 bends the auxiliary tube 58 and thereby shifts the main axis 76 of the surrounding expandable structure 56 relative to the axis 78 of the access path (i.e., the guide sheath 66). The prebent stylet 74, positioned within the interior of the structure 56, aids in altering the orientation of the structure 56 within targeted region. It is thereby possible to orient the structure 56 in a more generally aligned relationship with the natural axes of the vertebral body 26. A cavity 70, more centrally located within the bone, e.g., a vertebral body 26, can be established, which provides more uniform support across the mid region of the vertebral body 26 when filled with bone cement. The capability of the vertebral body 26 to withstand loads is thereby enhanced. The symmetric compaction of cancellous bone 32 in the interior volume also exerts more equal and uniform interior forces upon cortical bone 32, to elevate or push broken and compressed bone.

There are times when a lesser amount of cancellous bone compaction is indicated. For example, when the bone disease being treated is localized, such as in avascular necrosis, or where local loss of blood supply is killing bone in a limited area, an expandable structure 56 or 80 or 130 can compact a smaller volume of total bone. This is because the

diseased area requiring treatment is smaller.

Another exception lies in the use of an expandable structure 56 or 80 or 130 to improve insertion of solid materials in defined shapes, like hydroxyapatite and components in total joint replacement. In these cases, the structure shape and size is defined by the shape and size of the material being inserted.

Yet another exception lies in the use of expandable structures in bones to create cavities to aid in the delivery of therapeutic substances, as disclosed in copending United States Patent Application Serial No. 08/485,394, previously mentioned. In this case, the cancellous bone may or may not be diseased or adversely affected. Healthy cancellous bone can be sacrificed by significant compaction to improve the delivery of a drug or growth factor which has an important therapeutic purpose. In this application, the size of the expandable structure 56 or 80 or 130 is chosen by the desired amount of therapeutic substance sought to be delivered. In this case, the bone with the drug inside may need to be supported by standard methods while the drug works and the bone heals.

III. SINGLE USE

Distention of any one of the expandable structures 56 or 80 or 130 described herein during first use in a targeted body region generates stress on the material or materials which make up the structure. The material stress created by operational loads during first use in a targeted body region can significantly alter the preformed morphology of the structure, making future performance of the structure unpredictable.

For example, expansion within bone during a single use creates contact with surrounding cortical and cancellous bone. Regardless of the superior mechanical properties of material, this contact can in time damage the structure, creating localized regions of weakness, which may escape

5 detection. Localized areas of lower density cancellous bone may result in creating areas of differential expansion and stress on the structure. The existence of localized regions of weakness or differential stress can unpredictably cause overall structural failure during a subsequent use.

10 In addition, exposure to blood and tissue during a single use can entrap biological components on or within the structure or the associated catheter tube. Despite cleaning and subsequent sterilization, the presence of entrapped biological components can lead to unacceptable pyrogenic reactions.

15 As a result, following first use, the structure can not be consistently relied upon to reach its desired configuration during subsequent use and may not otherwise meet established performance and sterilization specifications. The effects of material stress and damage caused during a single use, coupled with the possibility of pyrogen reactions even after resterilization, reasonably justify imposing a single use restriction upon devices which carry these expandable structures for deployment in bone.

20 To protect patients from the potential adverse consequences occasioned by multiple use, which include disease transmission, or material stress and instability, or decreased or unpredictable performance, the invention also provides a kit 100 (see Figs. 14 and 15) for storing a single use tool 48 (also shown in Fig. 3) prior to use. As shown in Fig. 14, the tool 48 carries an expandable structure. Fig. 14 shows for the purpose of illustration the structure 56, as described herein. It should be appreciated that the tool 48 could carry an expandable structure 80 or 130, as also previously described.

30 In the illustrated embodiment (see Figs. 14 and 15), the kit 100 includes an interior tray 108. The tray 108 holds the tool 48 in a lay-flat, straightened condition during sterilization and storage prior to its first use.

The tray 108 can be formed from die cut cardboard or thermoformed plastic material. The tray 108 includes one or more spaced apart tabs 110, which hold the catheter tube 50 and expandable structure 56 in the desired lay-flat, straightened condition.

The kit 100 includes an inner wrap 112, which is peripherally sealed by heat or the like, to enclose the tray 108 from contact with the outside environment. One end of the inner wrap 112 includes a conventional peel-away seal 114, to provide quick access to the tray 108 upon instance of use, which preferably occurs in a sterile environment, such as within an operating room.

The kit 100 also includes an outer wrap 116, which is also peripherally sealed by heat or the like, to enclosed the inner wrap 112. One end of the outer wrap 116 includes a conventional peel-away seal 118, to provide access to the inner wrap 112, which can be removed from the outer wrap 116 in anticipation of imminent use of the probe 102, without compromising sterility of the probe 102 itself.

Both inner and outer wraps 112 and 116 (see Fig. 15) each includes a peripherally sealed top sheet 120 and bottom sheet 122. In the illustrated embodiment, the top sheet 120 is made of transparent plastic film, like polyethylene or MYLAR™ material, to allow visual identification of the contents of the kit 100. The bottom sheet 122 is made from a material that is permeable to ETO sterilization gas, e.g., TYVEK™ plastic material (available from DuPont).

The sterile kit 100 also carries a label or insert 106, which includes the statement "For Single Patient Use Only" (or comparable language) to affirmatively caution against reuse of the contents of the kit 100. The label 106 also preferably affirmatively instructs against resterilization of the tool 48. The label 106 also preferably instructs the physician or user to dispose of the tool 48 and the entire contents of the kit 100 upon use in accordance with

applicable biological waste procedures. The presence of the probe 102 packaged in the kit 100 verifies to the physician or user that tool 48 is sterile and has not been subjected to prior use. The physician or user is thereby assured that
5 the expandable structure 56 meets established performance and sterility specifications, and will have the desired configuration when expanded for use.

The label 106 preferably also instructs the physician as to the use of the expandable structure 56 (or 80 or 130)
10 for compacting cancellous bone in the manners previously described. For example, the label 106 instructs the physician to expand the structure inside bone to compact cancellous bone and form a cavity. The label 106 can also
15 instruct the physician to fill the cavity with a material, e.g., bone cement, allograft material, synthetic bone substitute, a medication, or a flowable material that sets to a hardened condition.

The features of the invention are set forth in the following claims.

09037250-041801